



Photo by Hadassa Koning

**“The development of a diagnostic tool in order to contain the
Vinegar Syndrome by effective collection management”**

Hadeel Abdul Aziz

IAQ Braunschweig, 15-16 November 2006



Content presentation

- Consortium
- Background
- Detection method
- Goal of the project
- Experimental
- Results
- Impressions
- Conclusions



Consortium

- ✓ TNO Indoor Health and Environment (lead contractor)
- ✓ Art Conservation bv
- ✓ Buchi Labortechnik GmbH
- ✓ Nederlands Filmmuseum
- ✓ Nationaal Archief



Background 1

Vinegar Syndrome is a threat!

- Employees, exposure to acetic acid → harmful to health
- Cultural heritage → deterioration photographic materials




Photo by Hadassa Koning

Background 2

Methods currently used to manage photographic collections:

- determination type of material: destructive or by visual observation
- observation of odour by nose (vinegar) in combination with AD-strips

Need for a rapid and objective method to determine:

- type of photographic material (cellulose acetate, cellulose nitrate or polyester)
- degradation state of the materials



Detection method

Detection method

1. Rapid
2. Reproducible
3. Objective

→ Method is based on Near Infrared spectroscopy

Method has successfully been applied in:

- Pharmaceuticals
- Food and nutrition
- Paper and pulp



Goal of the project

Application of Near Infrared Spectroscopy in order to determine the class and degradation state of photographic films.




Experimental

- Selection materials (consortium)
- Description materials
- Determination degradation state using AD-strips
- Recording NIR spectra
- Translation NIR spectra to models




Results 1

Model: Acetate / Nitrate / Polyester




Results 2



Results 3

Degradation model: Cellulose Acetate



Results 4

Example routine measurement

Input

Sample no.	123
Date	15/11/2006

Output


Sample 123	Cellulose acetate
Origin	DS 1-2
Degradation	




Impressions (1)



Impressions (2)



Impressions (3)



Conclusions

- The origin and degradation state of the films and transparencies can be easily performed using NIR
- Importance of this method.
 - Conventional methods were often failing to determine the origin of the photographic materials.
 - Even if senior conservators were using it.
 - This method proved to be a reliable and fast tool.

Acknowledgements

Maarten de Voer, Dr. Christoph Lühr and Dr.
Jochen Knecht Büchi Labortechnik

Piet van Dalen and Hadassa Koning Art
Conservation bv

Ad Pollé and Walter Swagemakers (Filmmuseum)

Gerrit de Bruin and Joop Korswagen (Nationaal
Archief)

Dr. John Havermans, Willem Suitela and Marco
Hubers (TNO)

